

SCF PURE CORTICO CANCELLOUS CRUNCH – STANDARD 1MM - 3.5MM



PRODUCED FROM

Trabecular and Cortical Bone

PRESERVATION/STORAGE METHOD

Frozen

COMMONLY USED FOR

Various Bone Void Applications, Spinal Surgery and Various Orthopaedic Procedures

- Over 15000 SCF Allografts implanted in Australia
- Consists of approx. 30% cortico and 70% cancellous bone granules
- Resorbs naturally, highly porous, highly permeable
- Ready to use in pure form or can be mixed with blood, bone marrow aspirate or saline
- Osteoconductive matrix promotes incorporation

DESCRIPTION	SPECIFICATIONS	CODE
SCF Cortico Cancellous Crunch	2.5g	AB-SCF2.5
SCF Cortico Cancellous Crunch	7.5g	AB-SCF7.5

- All donated tissue and allograft implants are screened and processed in accordance with the strict guidelines of the Australian Therapeutic Goods Administration
- SCF Allografts are processed and sterilized using Super Critical Fluid Technology
- All SCF Allografts are supplied as sterile (sterility assurance level 10^{-6}) using our patented release process