


TRICORTICAL CANCELLOUS WEDGE



PRODUCED FROM
Hemi Pelvis

PRESERVATION/STORAGE METHOD
Frozen

COMMONLY USED FOR
Anterior Cervical Fusion, Discectomy and
Various Orthopaedic Procedures

SPECIFICATIONS

CODE

Height	Width	Length	
10-20mm	10mm minimum	20-30mm	AB-TA108

*Variance on all measurements ± 1 mm

- All donated tissue and allograft implants are screened and processed in accordance with the strict guidelines of the Australian Therapeutic Goods Administration
- Allografts are manufactured utilizing our proprietary anti-bacterial rinses to ensure a final Bioburden of less than < 1 CFU
- All Allografts are treated to low dose gamma irradiation prior to release