


CORTICO- CANCELLOUS BLOCKS



PRODUCED FROM

Proximal Tibia, Proximal Femoral

PRESERVATION/STORAGE METHOD

Frozen

COMMONLY USED FOR

Acetabulum Repair, Total Hip/Knee Revision,
Maxillofacial Surgery and Other Procedures

SPECIFICATIONS

CODE

Width	Length	
15mm	60mm	AB-TA101

- *Cortical should be minimum 2mm thick
- *Cancellous bone maximum 3mm thick
- *Variance on all measurements = ± 1.0 mm

- All donated tissue and allograft implants are screened and processed in accordance with the strict guidelines of the Australian Therapeutic Goods Administration
- Allografts are manufactured utilizing our proprietary anti-bacterial rinses to ensure a final Bioburden of less than <1 CFU
- All Allografts are treated to low dose gamma irradiation prior to release