



# FEMORAL HEAD

**PRODUCED FROM**  
Proximal Femoral

**PRESERVATION/STORAGE METHOD**  
Frozen

**COMMONLY USED FOR**  
Total Joint Revisions,  
Various Orthopaedic Procedures

## SPECIFICATIONS

## CODE

Diameter

35mm minimum

AB-TA102

\*Variance on all measurements  $\pm 1$ mm

- All donated tissue and allograft implants are screened and processed in accordance with the strict guidelines of the Australian Therapeutic Goods Administration
- Allografts are manufactured utilizing our proprietary anti-bacterial rinses to ensure a final Bioburden of less than  $<1$  CFU
- All Allografts are treated to low dose gamma irradiation prior to release