



HTO WEDGE

PRODUCED FROM
Tibial Plateau

PRESERVATION/STORAGE METHOD
Frozen

COMMONLY USED FOR
High Tibial Osteotomy,
Various Orthopaedic Procedures

SPECIFICATIONS			CODE
Height	Width	Length	
10mm	Measure	50mm minimum	AB-TA104
12mm	Measure	50mm minimum	AB-TA105
14mm	Measure	50mm minimum	AB-TA106
16mm	Measure	50mm minimum	AB-TA107

*Measure: size variations in grafts are due to the nature of their source. Please contact Australian Biotechnologies for individual requirements.

*Variance on all measurements ± 1 mm

*Height is measured on highest side

- All donated tissue and allograft implants are screened and processed in accordance with the strict guidelines of the Australian Therapeutic Goods Administration
- Allografts are manufactured utilizing our proprietary anti-bacterial rinses to ensure a final Bioburden of less than <1 CFU
- All Allografts are treated to low dose gamma irradiation prior to release