




ACHILLES TENDON WITH BONE BLOCK



PRODUCED FROM

Achilles Tendon and Calcaneus Bone

PRESERVATION/STORAGE METHOD

Frozen

COMMONLY USED FOR

Anterior Cruciate Ligament Reconstruction,
Posterior Cruciate Ligament Reconstruction

SPECIFICATIONS

CODE

Graft	Width	Length	
Bone Block	12mm minimum	23mm minimum	AB-ST101
Tendon	—	180mm minimum	
Bundle	Measure	—	

*Variance on all measurements ± 1 mm

*Measure: size variations in grafts are due to the nature of their source. Please contact Australian Biotechnologies for individual requirements

- All donated tissue and allograft implants are screened and processed in accordance with the strict guidelines of the Australian Therapeutic Goods Administration
- Allografts are manufactured utilizing our proprietary anti-bacterial rinses to ensure a final Bioburden of less than <1 CFU
- All Allografts are treated to low dose gamma irradiation prior to release