




PATELLA TENDON WITH BONE BLOCK



PRODUCED FROM
Tibia, Patella

PRESERVATION/STORAGE METHOD
Frozen

COMMONLY USED FOR
Anterior Cruciate Ligament Reconstruction,
Posterior Cruciate Ligament Reconstruction

SPECIFICATIONS

CODE

Graft	Width	Length	
Bone Blocks	15mm minimum	25mm minimum	
Tendon	*same as width of bone block	Measure	AB-ST201
W / Quad	Measure	45mm minimum	Custom only

*Variance on all measurements ± 1 mm

*Measure: size variations in grafts are due to the nature of their source. Please contact Australian Biotechnologies for individual requirements

- All donated tissue and allograft implants are screened and processed in accordance with the strict guidelines of the Australian Therapeutic Goods Administration
- Allografts are manufactured utilizing our proprietary anti-bacterial rinses to ensure a final Bioburden of less than <1 CFU
- All Allografts are treated to low dose gamma irradiation prior to release