



# TIBIALIS TENDON

## ANTERIOR

**PRODUCED FROM**  
Anterior Tibialis Tendons

**PRESERVATION/STORAGE METHOD**  
Frozen

**COMMONLY USED FOR**  
Anterior/Posterior Cruciate Ligament  
Reconstruction

	SPECIFICATIONS		CODE
	Diameter	Length	
Anterior	7.5mm minimum (Folded)	230mm minimum	AB-ST301

\*Diameter of tendons is measured by folding tendon once and measuring combined diameter

\*Variance on all measurements  $\pm 1$ mm

- All donated tissue and allograft implants are screened and processed in accordance with the strict guidelines of the Australian Therapeutic Goods Administration
- Allografts are manufactured utilizing our proprietary anti-bacterial rinses to ensure a final Bioburden of less than  $<1$  CFU
- All Allografts are treated to low dose gamma irradiation prior to release