




TIBIALIS TENDON POSTERIOR



PRODUCED FROM
Posterior Tibialis Tendons

PRESERVATION/STORAGE METHOD
Frozen

COMMONLY USED FOR
Posterior Cruciate Ligament Reconstruction

	SPECIFICATIONS		CODE
	Diameter	Length	
Posterior	7.5mm minimum (Folded)	230mm minimum	AB-ST302

*Diameter of tendons is measured by folding tendon once and measuring combined diameter

*Variance on all measurements ± 1 mm

- All donated tissue and allograft implants are screened and processed in accordance with the strict guidelines of the Australian Therapeutic Goods Administration
- Allografts are manufactured utilizing our proprietary anti-bacterial rinses to ensure a final Bioburden of less than <1 CFU
- All Allografts are treated to low dose gamma irradiation prior to release