



# PERONEUS LONGUS TENDON

**PRODUCED FROM**  
Peroneus Longus Tendon

**PRESERVATION/STORAGE METHOD**  
Frozen

**COMMONLY USED FOR**  
Anterior Cruciate Ligament Reconstruction,  
Posterior Cruciate Ligament Reconstruction

SPECIFICATIONS		CODE
<b>Diameter</b>	<b>Length</b>	
7.5mm minimum (Folded)	230mm minimum	AB-ST401

\*Diameter of tendons is measured by folding tendon once and measuring combined diameter

\*Variance on all measurements  $\pm 1$ mm

- All donated tissue and allograft implants are screened and processed in accordance with the strict guidelines of the Australian Therapeutic Goods Administration
- Allografts are manufactured utilizing our proprietary anti-bacterial rinses to ensure a final Bioburden of less than  $<1$  CFU
- All Allografts are treated to low dose gamma irradiation prior to release