



### SCF PURE CORTICO CANCELLOUS CRUNCH – FINE 500 $\mu$ - 1MM

#### PRODUCED FROM

Trabecular and Cortical Bone

#### PRESERVATION/STORAGE METHOD

Frozen

#### COMMONLY USED FOR

Revision Bone Void Applications,  
Revision Orthopaedic Procedures and  
OMF Procedures

- Over 15000 SCF Allografts implanted in Australia
- Consists of approx. 30% cortico and 70% cancellous bone granules
- Resorbs naturally, highly porous, highly permeable
- Ready to use in pure form or can be mixed with blood, bone marrow aspirate or saline
- Osteoconductive matrix promotes incorporation

DESCRIPTION	SPECIFICATIONS	CODE
SCF Cortico Cancellous Crunch fine granules	2.5g	AB-SCFFG2.5
SCF Cortico Cancellous Crunch fine granules	1g	AB-SCFFG1

- All donated tissue and allograft implants are screened and processed in accordance with the strict guidelines of the Australian Therapeutic Goods Administration
- SCF Allografts are processed and sterilized using Super Critical Fluid Technology
- All SCF Allografts are supplied as sterile (sterility assurance level  $10^{-6}$ ) using our patented release process