



# CANCELLOUS CUBE

## PRODUCED FROM

Femoral Head, Femoral Condyle, Cancellous Bone

## PRESERVATION/STORAGE METHOD

Frozen

## COMMONLY USED FOR

Various Bone Void Applications,  
Various Dental and OMF procedures,  
Spinal Surgery

## SPECIFICATIONS

## CODE

Width	Length	Height	
10mm	10mm	10mm	AB-TA301

\*Variance on all measurements  $\pm 1$ mm

- All donated tissue and allograft implants are screened and processed in accordance with the strict guidelines of the Australian Therapeutic Goods Administration
- Allografts are manufactured utilizing our proprietary anti-bacterial rinses to ensure a final Bioburden of less than  $<1$  CFU
- All Allografts are treated to low dose gamma irradiation prior to release