

CANCELLOUS BLOCK



PRODUCED FROM

Femoral Head, Femoral Condyle

PRESERVATION/STORAGE METHOD

Frozen

COMMONLY USED FOR

Various Bone Void Applications,
Spinal Surgery

SPECIFICATIONS

CODE

Width	Length	Height	
10mm	20mm	10mm	AB-TA302

*Variance on all measurements ± 1 mm

- All donated tissue and allograft implants are screened and processed in accordance with the strict guidelines of the Australian Therapeutic Goods Administration
- Allografts are manufactured utilizing our proprietary anti-bacterial rinses to ensure a final Bioburden of less than < 1 CFU
- All Allografts are treated to low dose gamma irradiation prior to release