



# CANCELLOUS PLUG

**PRODUCED FROM**

Femoral Condyle, Femoral head, Proximal Tibia

**PRESERVATION/STORAGE METHOD**

Frozen

**COMMONLY USED FOR**

Various Bone Void Applications,  
Various Dental & OMF Procedures,  
Spinal Surgery

**SPECIFICATIONS**

**CODE**

Diameter	Length	
6mm	12mm	AB-TA303

\*Variance on all measurements = ±1mm

- All donated tissue and allograft implants are screened and processed in accordance with the strict guidelines of the Australian Therapeutic Goods Administration
- Allografts are manufactured utilizing our proprietary anti-bacterial rinses to ensure a final Bioburden of less than <1 CFU
- All Allografts are treated to low dose gamma irradiation prior to release