

HEMI FEMORAL HEAD



PRODUCED FROM
Proximal Femoral

PRESERVATION/STORAGE METHOD
Frozen

COMMONLY USED FOR
Joint Revisions,
Various Orthopaedic Procedures

SPECIFICATIONS

CODE

Diameter

35mm minimum

AB-TA103

*Variance on all measurements ± 1 mm

- All donated tissue and allograft implants are screened and processed in accordance with the strict guidelines of the Australian Therapeutic Goods Administration
- Allografts are manufactured utilizing our proprietary anti-bacterial rinses to ensure a final Bioburden of less than <1 CFU
- All Allografts are treated to low dose gamma irradiation prior to release