



PATELLA TENDON PRE-SHAPED



PRODUCED FROM
Tibia, Patella

PRESERVATION/STORAGE METHOD
Frozen

COMMONLY USED FOR
Anterior Cruciate Ligament Reconstruction,
Posterior Cruciate Ligament Reconstruction

SPECIFICATIONS

CODE

Bone Blocks

Tendon

Diameter

Length each

Width

Length

10mm

23mm min

minimum 10mm Measure

AB-ST203

*Variance on all measurements ± 1 mm

*Measure: size variations in grafts are due to the nature of their source. Please contact Australian Biotechnologies for individual requirements

- All donated tissue and allograft implants are screened and processed in accordance with the strict guidelines of the Australian Therapeutic Goods Administration
- Allografts are manufactured utilizing our proprietary anti-bacterial rinses to ensure a final Bioburden of less than <1 CFU
- All Allografts are treated to low dose gamma irradiation prior to release