



CERVICAL SPACER INSTRUMENTATION

DESCRIPTION

Customised Instrument set for Australian Biotechnologies' Cervical Allograft

USE

Anterior Cervical Fusion

SPECIFICATIONS

Allograft Insertion Tool	Allograft Delivery System
Tamp Tool	For Precise Placing of Allograft
Handle Piece	Used for Sizing Tools and Tamp Tool
Parallel Sizing Tools	5mm 6mm 7mm 8mm 9mm
Lordotic Sizing Tools	5mm 6mm 7mm 8mm 9mm

*Variance on all measurements ± 1 mm

*Height is measured on highest side